

MTS Personal Care Series

MTS ROLLER - CR02

1 Product Type & Info

Product Classification:	MTS Treatment Product
Product Type:	Microneedle Roller
Country of Origin:	U.S.A.
Available Size:	0.2mm



2 Description


It is a procedure which stimulates the skin to produce new collagen formation, thereby reducing wrinkles, and improving skin texture. A roller with very fine needles is used to produce many microscopic channels in the skin, which stimulate your own body to produce new collagen. These channels significantly improve the penetration of therapeutic active ingredients (called MTS Treatments) as well as multi-peptides, vitamins, or any other liposomal serums. MTS Treatments stimulate skin renewal, thereby making the skin appear fresher and younger. MTS is different to artificial injectable collagen, as the effects of collagen induction activities are long lasting.

MTS Roller is a durable and high-quality precision tool. The barrel and handle are made from Lexan® polycarbonate resin. The axle is polished stainless steel axle. The 200 bristles (needles) – the highest count of any roller – are made from surgical stainless steel. A patented manufacturing process holds the needles permanently in place to within 0.02mm of the specified length from the roller base. All rollers have unique quality-control serial numbers.

3 Key Points

- The treatment is almost free of pain
- There are no known negative side-effects
- You can use it at home, when ever you wish
- Apply your cosmetic on the skin – and roll it in
- A precision German-made instrument that will last for hundreds of treatments.
- The process stimulates physiological responses

5 Direction for use



Divide the treatment area(s) into sub-divisions, for instance the forehead, nose, cheeks, mouth area, and neck area. Apply MTS-Roller with the same amount of pressure on the skin throughout. Roll the same areas vertically and horizontally.

6 Cautions

Keep out of reach of children.