

# PRO ORDER PACKET

## Application Form and Professional Use Waiver

**MTS Roller™**

MTS Roller™ (Microneedle Therapy System™) is a skin treatment device designed and manufactured by - **Clinical Resolution Lab, Inc.** (a California corporation). All MTS Roller models are currently available for U.S. market (see below) and are FDA registered as Type I supplementary medical instruments. All components are produced in high-tech ISO13485 quality-controlled manufacturing facilities. MTS Rollers are made in U.S.A. and patented.

Pictured: Model MR20 (2.0mm)



**General Description**

MTS Roller is a durable and high-quality precision tool. The barrel and handle are made from medical grade Lexan® polycarbonate resin. The axle is polished stainless steel axle. The 200 bristles (needles) are made from surgical stainless steel. A patented manufacturing process holds the needles permanently in place to within 0.02mm of the specified length from the roller base. All rollers have unique quality-control serial numbers.

**GUARANTEE: 100% Contaminant-Free**

MTS Rollers are sterilized with gamma radiation and sealed in a hermetic pouch and are guaranteed to be 100% contaminate-free upon delivery.

**Storage and Use**

Polycarbonate resin does not degrade from prolonged use or storage. Roller models (MR2, MR3) should be replaced after ~6 months of regular use (2 to 3 times per week). A clinical roller (MR5) can be used up to ~10 treatments. Medical Roller (MR10, MR15, MR20) is only for one time use. Rollers that have been dropped or bumped should be replaced immediately. MTS Rollers should not be transferred between patients. Keep rollers in the protective case when not in use. For more information, consult CRL protocol guides.

**MTS CLINICAL SOLUTIONS**  
(Helps to Enhance MTS Roller Treatment Results)



**Care of MTS Rollers™**

Current MTS Roller models can withstand repeated immersion in “cold sterilization” agents such as germicides (phenol, ethyl alcohol) and Cidex™ (glutaraldehyde). Rollers should be completely dry before soaking to prevent dilution of the sterilization agent. Soak time depends upon the manufacturer's specifications. Rinse rollers thoroughly before using and return to its protective case after drying.



<b>MTS ROLLER MODELS:</b>	
<b>MR2:</b>	<b>0.2mm (personal use)</b>
<b>MR3:</b>	<b>0.3mm (personal/ clinical use)</b>
<b>MR5:</b>	<b>0.5mm (clinical use only)</b>
<b>MR10:</b>	<b>1.0mm (medical use only)</b>
<b>MR15:</b>	<b>1.5mm (medical use only)</b>
<b>MR20:</b>	<b>2.0mm (medical use only)</b>

# PRO ORDER PACKET

Application Form and Professional Use Waiver

## To Order MTS Roller™

Complete the below information contained in this “Pro Order Packet” as a convenient means to begin a unique opportunity for professional clinicians (physicians and/or other registered/licensed/certified skincare professionals) to test MTS Roller in their practice at DEALER PRICES. Please FAX BACK completed form along with a copy of your PROFESSIONAL LICENSE to our office - **213-365-0053** – or email a PDF/Scanned version of the completed document to [info@clinicalresolution.com](mailto:info@clinicalresolution.com). Please note that there are TWO parts to this Packet: “Certificate of Professional Clinical Status,” and “Waiver of Liability”. For any questions please call (877) KNOW-MTS (566-9687).

### Certification of Professional Clinical Status

Practitioner’s Name (Buyer):	Phone:
Business Name (dba):	Fax:
Mailing/Billing Address (stamp or label):	Email/Website
	License/Certification Number: <i><b>*Please fax a copy of your License/Certificate along with this Pro Order Packet</b></i>
	Are you currently in practice?    YES or NO
Shipping Address (if different than mailing/billing):	Years in practice
	Your Clinic Specialty(ies)
Reseller Permit # (if in California)	I am the “Practitioner” named above. All of the information I have provided is accurate to the best of my ability. (Sign and date)

### Waiver of Liability

Buyer (named signor on the “Certificate of Professional Clinical Status” and also signing below) assumes sole responsibility for the use of MTS Rollers. Buyer also assumes complete and full responsibility – without exception – for any liability (material, financial, medical, and/or psychological) that may arise from the use of MTS Rollers and/or solutions. Buyer completely indemnifies and holds harmless Clinical Resolution Lab, Inc., its officers, representatives, and/or agents from any judgment, any arbitration, and/or any other settlement brought with or without cause.

Clinical Resolution Lab, Inc. will replace defective MTS Rollers that fail during normal usage. Buyer also acknowledges that until further notice, Clinical Resolution advises that MTS Rollers MR2, MR3 and MR5 should be sterilized by cold sterilization processes between usages. **Do not autoclave MTS Rollers.** This warranty is therefore exempt for MTS Rollers that have been autoclaved.

MTS Rollers MR10, MR15 and MR20 are sold as single-patient use supplemental surgical tools for use in the clinical environment only under supervision of a physician. MTS Rollers MR2 and MR3 are sold as single-patient use take-home tools to be used under the guidance of a physician.

Other terms and conditions apply per CRL, Inc. invoicing and payment policies.

*I agree to the terms and conditions of the above Waiver of Liability. (sign and date)*

Signature:	Date:
------------	-------