MTS Personal Care Series

MTS HAIR RESTORATION COMPLEX

Product Type & Info

Product Classification: Product Type: Manufactured by: Available Size:

MTS Hair Regrowth Solution Package MTS 0.2mm & Hair product Clinical Resolution Lab., Inc. 1 Box



2 Description

The most effective MTS process for transdermal delivery of cosmetics into deeper skin layers.

In contrast to normal cosmetics, MTS Anti-Aging Complex offers a totally innovative way to rejuvenating your skin. The MTS-Roller (0.2mm version inside) enhances our specially designed MTS Anti-Aging Complex cosmetics into deeper skin layers - where all the skin's problems start.

MTS Peptide HR Complex is the first and only product in its class to directly stimulate the Anagenic Growth phase of hair cycle with growth promoting growth factors for hair.

The application of these findings through our advances in peptide biochemistry offers a new beginning in hair restoration technology for both men and women.

3 Contents

1	MTS Peptide HR Matrix
1	MTS Roller

(CR2: 0.2mm)

(3.3 fl.oz/100ml) Poly-Vitamin and Peptide Treatment New Era of Skin Rejuvenation with 200 microneedles with GE Lexan polycarbonate resin

- MTS Roller Cleaning Tray 1
- MTS Roller Cleaning Brush 1
- **MTS Sterilizing Solution** 1 (6 fl.oz/180ml) Isopropyl Alcohol 70%, Water

4 Key Points

MTS Hair Restoration Complex improves:

- At 8 weeks, 100% of men treated with MTS Peptide HR Matrix experienced an increase in hair growth and an average 80% reduction in hair loss.
- Twelve weeks of study revealed a 62% reduction in sebum production for volunteers using MTS Peptide HR Complex.

5 Direction for use

- 1. Step 1 Boosting: MTS VitaPep Serum
 - MTS-Roller
- 2. Step 2 Correction: 3. Step 3 - Maximum Effect: Repeat Step 1
- 4. Step 4 Cleaning: MTS Sterilizing Solution.

